

SPRI Bead Extraction:
Version 6 (13 July 2020)

Notes:

- This protocol has been optimized using nucleated avian blood stored on lysis buffer.
- Beads **MUST** be at room temp and well mixed
 - Set out from fridge ~30 mins prior to starting
 - Also a good time to set out TLE (elution buffer) if that's what you'll use
- You must prepare **FRESH** 80% EtOH. This is critical because if too low and the DNA will elute off the beads, and too high you are likely to over-dry due to evaporation.
- Set the incubator to 37C while you're getting set up.

Digestion:

1. Use wide-bore 1,000uL pipette to wrangle ~200-400uL starting material into a new 1.5uL Eppendorf
 - a. Be sure you're using low-bind tubes!
 - b. May have to "cut" the blood blob by closing cap
2. Add 30uL of 20 mg/mL Prot-K and 20uL 10X Queen's Lysis buffer
 - a. Vortex for 5 seconds to mix and flash spin to collect any drops
3. Incubate at 66C / 250 RPM overnight in the shaking incubator

DNA Isolation:

1. Combine **200uL of digest with 300uL SpeedBeads**
 - a. **MUST** BE 1-2X beads than digest.
 - b. Save extra digest product in case extraction fails
 - c. Mix by pipetting with **WIDE BORE** tips, spin down only if necessary
 - d. Let tubes sit for **5 mins. at room temperature**
2. Place **tubes on magnet for 15-30 mins** or until the solution is clear
 - a. Be patient! All your DNA is binding to the magnet right now so give it time!
 - b. **While you wait, make the 80% EtOH and TLE aliquot.**
3. **Remove and dispose supernatant** without disturbing the beads.
 - a. You **CAN KEEP** the supernatant at this point and do another round if needed
4. **Wash 1:** Add 700uL of **FRESH** 80% EtOH to the tube
 - a. This volume must be greater than volume of beads + digest solution
 - b. Let it sit for 1 min
 - i. Use a stopwatch if it's only a couple samples, but if you have ~10+ it will be over a minute before you get back around
 - c. Remove and dispose EtOH without disturbing the beads.

5. **Wash 2:** Add 800uL of 80% EtOH to the tube
 - a. Let it sit for 1 min.
 - b. Remove and dispose EtOH without disturbing the beads.
 - c. Take tubes off the magnet and flash spin (5 secs) to pellet beads
 - d. Place back on magnet and use 200uL pipette to remove remaining EtOH
6. **Add 200uL of TLE buffer** to each tube
 - a. Mix by flicking or gentle vortexing, flash spin if necessary
 - i. Need beads and DNA to disassociate, but don't want to shear DNA
 - ii. I found higher volume of TLE helps beads disassociate
7. **Incubate at 37C for ~15 mins.**
 - a. Keep tubes open so any EtOH carryover from the wash steps can evaporate
 - b. Use this time to make qubit buffer if that's next
8. Place **tubes back on magnet for ~10 mins** or until solution is clear
 - a. VERY IMPORTANT! The solution MUST be clear. If the DNA solution is highly concentrated this could take a while... Be patient!
9. **Remove supernatant (aka liquid gold)** and place into a new, labelled tube!
 - a. Cross fingers that everything worked ;)